

# FLUXANA

Your supplier for  
XRF Application Solutions

# HD ELEKTRONIK

Fusion machines

## Liquid Analysis Technique

We deliver everything you need for a professional liquid analysis.



### ■ Our new professional sample cup series

| Part No.            | Description   | OD mm | ID mm | H mm | Vol. ml |
|---------------------|---|-------|-------|------|---------|
| SC-3332             | Cup, ring, cap for: Spectro XEPOS, XLAB 2000®                     | 32    | 25    | 23   | 10      |
| SC-3335             | Cup, ring, cap for Philips/Panalytical®                           | 35    | 28    | 30   | 18      |
| SC-3340             | Inner ring, outer ring, cap for: Bruker®SRS, S4, S8, ARL Optimix® | 40    | 34    | 21   | 18      |
| <b>NEW</b> SC-3340T | Inner ring, outer ring, cap for: Bruker®SRS, S4, S8               | 40    | 34    | 29   | 25      |
| SC-3345             | Inner ring, outer ring, cap universal for 45mm                    | 45    | 37    | 40   | 40      |
| SC-3351             | Cup, cap universal for 51,5mm holder                              | 51    | 44    | 33   | 50      |

Please ask for a quote.



**FLUXANA GmbH & Co. KG**  
Sommerdeich 22, D-47533 Kleve, Germany  
Tel.: +49 (0) 2821 97 38 75  
Fax: +49 (0) 2821 97 38 76  
E-mail: info@fluxana.de  
Web: www.fluxana.com

Amtsgericht Kleve: HR-A 2935, HR-B 8211  
Ust-IdNr.: DE 814692564, Steuer-Nr. 116/5755/0442  
Finanzamt Kleve

**In cooperation with:**

**HD Elektronik und Elektrotechnik GmbH**  
Sommerdeich 22, D-47533 Kleve, Germany  
Tel.: +49 (0) 2821 148 10  
Fax: +49 (0) 2821 148 09  
E-mail: hde@hdelektronik.de  
Web: www.hdelektronik.de

Amtsgericht Kleve: HR-B 2162  
Ust-IdNr.: DE 812941185, Steuer-Nr. 116/5707/1869  
Finanzamt Kleve

Official agent

## Professional Cups

The sample cups of FLUXANA are designed for use in XRF. All cups have been optimized and approved by our customers world wide.



We offer different XRF-sample cups to you.

| Part No.     | Description  | OD mm | ID mm | H mm | Vol ml |
|--------------|--|-------|-------|------|--------|
| SC-4031      | Cup, ring, single open ended, universal open for 32mm holder   | 32    | 25    | 24   | 12     |
| SC-4047      | Cup, ring, single open ended, for Horiba®                      | 47    | 38    | 18,5 | 18     |
| SC-4131      | Cup, ring, collar single open ended, universal for 32mm holder | 32    | 25    | 24   | 12     |
| SC-4140      | Cup, ring, collar single open ended, universal for 40mm holder | 40    | 32    | 24   | 18     |
| SC-4231      | Cup, 2 rings, double open ended, universal for 32mm holder     | 32    | 25    | 24   | 12     |
| SC-4240      | Cup, 2 rings, double open ended, universal for 40mm holder     | 40    | 32    | 24   | 18     |
| SC-4331      | Cup, ring, cap universal for 32mm holder                       | 32    | 25    | 24   | 12     |
| SC-4335      | Cup, ring, cap for Philips/Panalytical®                        | 35    | 27    | 35   | 18     |
| SC-4340      | Cup, ring, cap universal for 40mm holder                       | 40    | 32    | 24   | 18     |
| SC-4345      | Cup, cap for Philips/Panalytical®                              | 45    | 39    | 32   | 36     |
| SC-6232      | Cup, 2 rings, double open ended for Spectro XEPOS, XLAB 2000®  | 32    | 25    | 23   | 10     |
| SC-6332      | Cup, ring, cap for Spectro XEPOS, XLAB 2000®                   | 32    | 25    | 23   | 10     |
| SC-7332      | Inner cup and cap, designed for Oxford®                        | 32    | 30    | 38,5 | 25     |
| SC-7332inner | Only inner cup, designed for Oxford LABX®                      | 32    | 30    | 38,5 | 25     |
| SC-7332cap   | Only cap, designed for Oxford LABX®                            |       |       |      |        |
| SC-8047      | Cup, ring, single open ended, for Horiba® 47mm holder          | 47    | 38    | 18,5 | 18     |

One set contains 100 pieces. Caps and rings are packed separately! No cutting necessary!



## Thin Films

The preferences of our thin films:

- 25% more pre-cut films per package than other brands
- dust proof anti static packaging
- Lot number tracking system
- rolls, pre-cut circles & sheets



| Part No.      | Description                             | Thickness µm       | Packaging                       |
|---------------|---|--------------------|---------------------------------|
| TF-025        | Microporous                             |                    | Roll 6,4cm x 5,1m               |
| TF-112        | Mylar                                   | 12                 | Roll 7,6cm x 91,4m              |
| TF-115        | Ultrapolyester                          | 1,5                | Roll 7,6cm x 91,4m              |
| TF-125        | Mylar                                   | 2,5                | Roll 7,6cm x 91,4m              |
| TF-125-345    | Mylar                                   | 2,5                | Square 7,6cm x 7,6cm (500 pcs.) |
| TF-12560      | PolyRep (replacement for polycarbonate) | Mylar 2,5 + PP 6,0 | 2 Rolls a 7,6cm x 91,4m         |
| TF-1335       | Hostaphan                               | 3,5                | Roll 7,6cm x 91,4m              |
| TF-1335-345   | Hostaphan                               | 3,5                | Square 7,6cm x 7,6cm (500 pcs.) |
| TF-135        | Mylar                                   | 3,5                | Roll 7,6cm x 91,4m              |
| TF-135-255    | Mylar                                   | 3,5                | Circle 6,4cm (500 pcs.)         |
| TF-135-345    | Mylar                                   | 6,0                | Square 7,6cm x 7,6cm (500 pcs.) |
| TF-160        | Mylar                                   | 6,0                | Roll 7,6cm x 91,4m              |
| TF-160-255    | Mylar                                   | 6,0                | Circle 6,4cm (500 pcs.)         |
| TF-160-345    | Mylar                                   | 6,0                | Square 7,6cm x 7,6cm (500 pcs.) |
| TF-212        | Polypropylene                           | 12,0               | Roll 7,6cm x 91,4m              |
| TF-212-345    | Polypropylene                           | 12,0               | Square 7,6cm x 7,6cm (500 pcs.) |
| TF-240        | Polypropylene                           | 4,0                | Roll 7,6cm x 91,4m              |
| TF-240-255    | Polypropylene                           | 4,0                | Circle 6,4cm (500 pcs.)         |
| TF-240-30     | Polypropylene                           | 4,0                | Circle 7,6cm (100 pcs.)         |
| TF-240-345    | Polypropylene                           | 4,0                | Square 7,6cm x 7,6cm (500 pcs.) |
| TF-250        | Polypropylene                           | 5,0                | Roll 7,6cm x 91,4m              |
| TF-260        | Polypropylene                           | 6,0                | Roll 7,6cm x 91,4m              |
| TF-260-255    | Polypropylene                           | 6,0                | Circle 6,4cm (500 pcs.)         |
| TF-260-345    | Polypropylene                           | 6,0                | Square 7,6cm x 7,6cm (500 pcs.) |
| TF-350        | Polycarbonate                           | 5,0                | Roll 7,6cm x 91,4m              |
| TF-475        | Kapton                                  | 7,5                | Roll 7,6cm x 91,4m              |
| TF-475-345    | Kapton                                  | 7,5                | Square 7,6cm x 7,6cm (500 pcs.) |
| TF-500        | XRF Tape                                |                    | Roll 25mm x 66m                 |
| TF-LS-240-255 | Polypropylene Low Sulfur                | 4,0                | Circle 6,4cm (500 pcs.)         |

